

Mobility gallery – Long-term care



A



Albert

- Ambulatory, but may use a walking stick for support
- Independent – can clean and dress himself
- Usually no risk of dynamic or static overload for the caregiver
- Stimulation of functional mobility is very important to maintain independence



B



Barbara

- Can support herself and may use a walking frame or similar
- Dependent on the caregiver in some situations
- Usually no risk of dynamic overload for the caregiver
- A risk of static overload may occur for example during assistance with daily hygiene activities
- Stimulation of functional mobility is very important



C



Carl

- Is able to partially weight bear on at least one leg. Often sits in a wheelchair and has some trunk stability
- Dependent on caregiver in most situations
- A risk of dynamic and static overload to the caregiver when not using proper equipment
- Stimulation of functional mobility for Carl



D



Doris

- Cannot stand and is not able to weight bear through her feet. Is able to sit if well supported
- Dependent on caregiver in most situations
- A high risk of dynamic and static overload to the caregiver when not using proper equipment
- Stimulation of functional mobility is very important



E



Emma

- Might be almost completely bedridden, can sit out only in a special chair
- Always dependent on caregiver
- A high risk of dynamic and static overload to the caregiver when not using proper equipment
- Stimulation of functional mobility is not a primary goal

Mobility gallery – Bariatric care



A



Albert

- Ambulatory, but may use a walking stick for support
- Independent – can clean and dress himself
- Usually no risk of dynamic or static overload for the caregiver
- Stimulation of functional mobility is very important to maintain independence



B



Barbara

- Can support herself and may use a walking frame or similar
- Dependent on the caregiver in some situations
- Usually no risk of dynamic overload for the caregiver
- A risk of static overload may occur for example during assistance with daily hygiene activities
- Stimulation of functional mobility is very important



C



Carl

- Is able to partially weight bear on at least one leg. Often sits in a wheelchair and has some trunk stability
- Dependent on caregiver in most situations
- A risk of dynamic and static overload to the caregiver when not using proper equipment
- Stimulation of functional mobility for Carl



D



Doris

- Cannot stand and is not able to weight bear through her feet. Is able to sit if well supported
- Dependent on caregiver in most situations
- A high risk of dynamic and static overload to the caregiver when not using proper equipment
- Stimulation of functional mobility is very important



E



Emma

- Might be almost completely bedridden, can sit out only in a special chair
- Always dependent on caregiver
- A high risk of dynamic and static overload to the caregiver when not using proper equipment
- Stimulation of functional mobility is not a primary goal

Therapeutic Seating

Think - 'Patients who use Sara Steady'

Criteria - For use with patients with limited mobility.

- Patients may require assistance with transfers
- Patient rehabilitation
- Risk of pressure injuries



Milano

- Weight limit: 350 lbs
- From 5° forward to 35° backward Tilt in Space
- Independent leg elevation and back recline
- Easy seat width and seat depth adjustment from 16" to 20"
- 3 Tier waterfall back or lateral back option
- Height adjustable footplate
- Fully manual and fully motorized options
- Sara Steady compatible
- Removable cushion



Think - 'The Falls Prevention Chair'

Criteria - For use with patients at high risk of sliding and falls from the chair.

- Patients with vigorous involuntary movements
e.g. Huntington's Disease
- Patients with Dementia
- Patients displaying agitation
- May require assistance with transfers



Atlanta

- Weight limit: 350 lbs
- Integrated tilt and recline for deep seat position
- Retractable transfer footplate
- Variable angle leg rest
- Compatible with patient lifters
- Removable cushion



Think - 'Postural support for those at high risk of pressure injuries'

Criteria - For use with patients with postural and seating needs who require upper body and head support.

- Patients with neurological conditions
e.g., Stroke, acquired brain injury
- Require full assistance with transfers
- Require upper body and head support
- Patient rehabilitation
- High risk of pressure injuries



Sorrento

- Weight limit: 350 lbs
- From 5° forward to 35° or 45° backward Tilt in Space
- Independent leg elevation and back recline
- Height adjustable footplate
- Removable arm for side transfers
- Lateral supports available
- Compatible with patient lifters
- Removable cushion

